



SpeciVar® BackLit Trays

- ▼ Speed Specimen Collection & Visualization
- ▼ Improve Ergonomics & Staff Satisfaction



Innovating GI Endoscopy takes more than a bright idea

Our vision of a successful solution doesn't stop at successful examinations. It starts with continually innovating products and improving staff work environments. The SpeciVar® BackLit Tray (BLT) is uniquely designed to bring this vision to light.

Speed Specimen Collection & Visualization

We recognize that specimen collection and visualization is a critical part of successful patient outcomes. Using a SpeciVar BLT speeds specimen collection and visualization by providing gentle back light in dimly lit procedure rooms.

Improve Ergonomics & Staff Satisfaction

The SpeciVar BLT's unique design* offers several staff benefits:


- Reduced eye strain that results from specimen collection in dim lighting
- A small footprint and cordless design preserves cart field space
- No additional turnover time is required since a BLT is easily wiped clean or can be fully submerged if needed
- A single or multiple well linear array prevents mislabeling specimen jars
- A low center of gravity creates a stable base reducing the risk of spills vs. freestanding specimen jars
- LED light sources remain cool after hours of use
- The light source intensity can be adjusted according to staff preference
- The light source can be changed in seconds with no special tools

There are three SpeciVar BLT models

Each backlit specimen tray:

- includes two USB rechargeable LED light sources with a total combined run time of 20 hrs
- is constructed using welded acrylic making it lightweight yet durable
- is designed to preserve cart field space and prevent spills

SPECIVUE
ENDOSCOPY ILLUMINATED



*SpeciVar® BackLit Trays are a protected patent pending design


To order or for more information visit:
www.SpeciVar.com

Staff will more than see the difference SpeciVar® BackLit Trays create in speed and improved ergonomics


With the SpeciVar® BackLit Tray (BLT)



Without the SpeciVar® BackLit Tray (BLT)



Minimalist, cordless design preserves cart field space



A multiple well linear array keeps specimen jars in order preventing mislabeling errors



LED backlit specimen jars make collection and visualization faster and easier reducing eye strain



Each BLT easily wipes clean between cases reducing turnover time



Wells are cut to match your specimen jars ensuring an easy fit into your existing process



Low center of gravity creates a stable base reducing spillage risks



To order or for more information visit:
www.SpeciVar.com



www.SpeciVue.com

